

# Introductory Lecture Series: The Anesthesia Machine

PORNSIRI WANNADILOK




## Objectives

- Anesthesia Machine
- Ventilators
- Scavenging Systems
- System Checkout



# The Anesthesia Machine



# High Pressure System



- Receives gasses from the high pressure E cylinders attached to the back of the anesthesia machine (2200 psig for O<sub>2</sub>, 745 psig for N<sub>2</sub>O)
- Consists of:
  - Hanger Yolk (reserve gas cylinder holder)
  - Check valve (prevent reverse flow of gas)
  - Cylinder Pressure Indicator (Gauge)
  - Pressure Reducing Device (Regulator)
- Usually not used, unless pipeline gas supply is off

## E Size Compressed Gas Cylinders

Cylinder Characteristics	Oxygen	Nitrous Oxide	Carbon Dioxide	Air
Color	White (green)	Blue	Gray	Black/White (yellow)
State	Gas	Liquid and gas	Liquid and gas	Gas
Contents (L)	625	1590	1590	625
Empty Weight (kg)	5.90	5.90	5.90	5.90
Full Weight (kg)	6.76	8.80	8.90	
Pressure Full (psig)	2000	750	838	1800

## Hanger Yolk

- Hanger Yolk: orients and supports the cylinder, providing a gas-tight seal and ensuring a unidirectional gas flow into the machine
- Index pins: Pin Index Safety System (PISS) is gas specific → prevents accidental rearrangement of cylinders (e.g.. switching O<sub>2</sub> and N<sub>2</sub>O)



## Pressure Reducing Device

- Reduces the high and variable pressures found in a cylinder to a lower and more constant pressure found in the anesthesia machine (45 psig)
- Reducing devices are preset so that the machine uses only gas from the pipeline (wall gas), when the pipeline inlet pressure is 50 psig.




This prevents gas use from the cylinder even if the cylinder is left open (i.e. saves the cylinder for backup if the wall gas pipeline fails)

## Pressure Reducing Device

- Cylinders should be kept closed routinely. Otherwise, if the wall gas fails, the machine will automatically switch to the cylinder supply without the anesthetist being aware that the wall supply has failed (until the cylinder is empty too).



## Intermediate Pressure System

- Receives gasses from the regulator or the hospital pipeline at pressures of 40-55 psig
- Consists of:
  - Pipeline inlet connections
  - Pipeline pressure indicators
  - Piping
  - Gas power outlet
  - Master switch
  - Oxygen pressure failure devices
  - Oxygen flush
  - Additional reducing devices
  - Flow control valves



## Pipeline Inlet Connections

- Mandatory N2O and O2, usually have air and suction too
- Inlets are non-interchangeable due to specific threading as per the Diameter Index Safety System (DISS)
- Each inlet must contain a check valve to prevent reverse flow (similar to the cylinder yolk)



## Oxygen Pressure Failure Devices


- Machine standard requires that an anesthesia machine be designed so that whenever the oxygen supply pressure is reduced below normal, the oxygen concentration at the common gas outlet does not fall below 19%

## Oxygen Pressure Failure Devices

- A **Fail-Safe valve** is present in the gas line supplying each of the flowmeters except O<sub>2</sub>. This valve is controlled by the O<sub>2</sub> supply pressure and shuts off or proportionately decreases the supply pressure of all other gasses as the O<sub>2</sub> supply pressure decreases
- Historically there are 2 kinds of fail-safe valves
  - Pressure sensor shut-off valve (Ohmeda)
  - Oxygen failure protection device (Drager)


## Pressure Sensor Shut-Off Valve

- Oxygen supply pressure opens the valve as long as it is above a pre-set minimum value (e.g.. 20 psig).
- If the oxygen supply pressure falls below the threshold value the valve closes and the gas in that limb (e.g.. N<sub>2</sub>O), does not advance to its flow-control valve.



## Oxygen Failure Protection Device (OFPD)

- Based on a proportioning principle rather than a shut-off principle
- The pressure of all gases controlled by the OFPD will decrease proportionately with the oxygen pressure



## Oxygen Supply Failure Alarm

- The machine standard specifies that whenever the oxygen supply pressure falls below a manufacturer-specified threshold (usually 30 psig) a medium priority alarm shall blow within 5 seconds.



## Limitations of Fail-Safe Devices/Alarms


- Fail-safe valves **do not** prevent administration of a hypoxic mixture because they depend on pressure and not flow.
- These devices **do not** prevent hypoxia from accidents such as pipeline crossovers or a cylinder containing the wrong gas

## Limitations of Fail-Safe Devices/Alarms

- These devices **prevent hypoxia from some problems occurring upstream** in the machine circuitry (disconnected oxygen hose, low oxygen pressure in the pipeline and depletion of the oxygen cylinder)
- Equipment **problems that occur downstream** (for example leaks or partial closure of the oxygen flow control valve) **are not prevented** by these devices.

## Oxygen Flush Valve (O<sub>2</sub> +)

- Receives O<sub>2</sub> from pipeline inlet or cylinder reducing device and directs high, unmetered flow directly to the common gas outlet (downstream of the vaporizer)
- Machine standard requires that the flow be between 35 and 75 L/min
- The ability to provide jet ventilation
- Hazards
  - May cause barotrauma
  - Dilution of inhaled anesthetic



## Second-Stage Reducing Device


- Located just upstream of the flow control valves
- Receives gas from the pipeline inlet or the cylinder reducing device and reduces it further to 26 psig for N<sub>2</sub>O and 14 psig for O<sub>2</sub>
- Purpose is to eliminate fluctuations in pressure supplied to the flow indicators caused by fluctuations in pipeline pressure

# Low Pressure System

- Extends from the flow control valves to the common gas outlet
- Consists of:
  - Flow meters
  - Vaporizer mounting device
  - Check valve
  - Common gas outlet


## Flowmeter assembly

- When the flow control valve is opened the gas enters at the bottom and flows up the tube elevating the indicator
- The indicator floats freely at a point where the downward force on it (gravity) equals the upward force caused by gas molecules hitting the bottom of the float




# Arrangement of the Flow-Indicator Tubes

- In the presence of a flowmeter leak (either at the "O" ring or the glass of the flow tube) a hypoxic mixture is less likely to occur if the O<sub>2</sub> flowmeter is downstream of all other flowmeters
- In A and B a hypoxic mixture can result because a substantial portion of oxygen flow passes through the leak, and all nitrous oxide is directed to the common gas outlet
- \* Note that a leak in the oxygen flowmeter tube can cause a hypoxic mixture, even when oxygen is located in the downstream position



# Proportioning Systems

- Mechanical integration of the N<sub>2</sub>O and O<sub>2</sub> flow-control valves
- Automatically intercedes to maintain a minimum 25% concentration of oxygen with a maximum N<sub>2</sub>O:O<sub>2</sub> ratio of 3:1



## Limitations of Proportioning Systems

- Machines equipped with proportioning systems can still deliver a hypoxic mixture under the following conditions:
  - Wrong supply gas
  - Defective pneumatics or mechanics (e.g.. The Link-25 depends on a properly functioning second stage regulator)
  - Leak downstream (e.g.. Broken oxygen flow tube)
  - Inert gas administration: Proportioning systems generally link only N<sub>2</sub>O and O<sub>2</sub>

## Vaporizers

- A vaporizer is an instrument designed to change a liquid anesthetic agent into its vapor and add a controlled amount of this vapor to the fresh gas flow



# Classification of Vaporizers

Methods of regulating output concentration

Concentration calibrated (e.g. variable bypass)

Measured flow

Method of vaporization

Flow-over

Bubble through

Injection

Temperature compensation

Thermocompensation


Supplied heat

## Generic Bypass Vaporizer

- Flow from the flowmeters enters the inlet of the vaporizer

- The function of the concentration control valve is to regulate the amount of flow through the bypass and vaporizing chambers

Splitting Ratio = flow through vaporizing chamber / flow through bypass chamber



## Factors That Influence Vaporizer Output


- **Flow Rate:** The output of the vaporizer is generally less than the dial setting at very low (< 200 ml/min) or very high (> 15 L/min) flows
- **Temperature:** Automatic temperature compensating mechanisms in bypass chambers maintain a constant vaporizer output with varying temperatures
- **Back Pressure:** Intermittent back pressure (e.g. positive pressure ventilation causes a higher vaporizer output than the dial setting)

## Factors That Influence Vaporizer Output

- **Atmospheric Pressure:** Changes in atmospheric pressure affect variable bypass vaporizer output as measured by volume % concentration, but not (or very little) as measured by partial pressure (lowering atmospheric pressure increases volume % concentration and vice versa)
- **Carrier Gas:** Vaporizers are calibrated for 100% oxygen. Carrier gases other than this result in decreased vaporizer output.

# The Circuit: Circle System

- Arrangement is variable, but to prevent re-breathing of CO<sub>2</sub>, the following rules must be followed:
  - Unidirectional valves between the patient and the reservoir bag
  - Fresh-gas-flow cannot enter the circuit between the expiratory valve and the patient
  - Adjustable pressure-limiting valve (APL) cannot be located between the patient and the inspiratory valve



## Circle System

- **Advantages:**
  - Relative **stability** of inspired concentration
  - Conservation of respiratory **moisture and heat**
  - Prevention of operating room **pollution**
  - PaCO<sub>2</sub> depends only on ventilation, not fresh gas flow
  - Low fresh gas flows can be used
- **Disadvantages:**
  - **Complex** design = potential for malfunction
  - High **resistance** (multiple one-way valves) = higher work of breathing





# The Adjustable Pressure Limiting (APL) Valve

- User adjustable valve that releases gases to the scavenging system and is intended to provide control of the pressure in the breathing system
- **Bag-mask Ventilation:** Valve is usually left partially open. During inspiration the bag is squeezed pushing gas into the inspiratory limb until the pressure relief is reached, opening the APL valve.
- **Mechanical Ventilation:** The APL valve is excluded from the circuit when the selector switch is changed from manual to automatic ventilation





# Scavenging Systems

- Protects the breathing circuit or ventilator from excessive positive or negative pressure.



# Scavenging Systems



## Checking Anesthesia Machines

8 Categories of check:

- Emergency ventilation equipment
- High-Pressure system
- Low-Pressure system
- Scavenging system
- Breathing system
- Manual and automatic ventilation system
- Monitors
- Final Position

