

DNA CLONING

Dr.Sarookhani





FIGURE 18.39 An example of DNA cloning using bacterial plasmids.

Dr.Sarookhani



Dr.Sarookhani




Fig. 6.1 Generalized overview of cloning strategies, with favored routes shown by arrows. Note that in cell-based cloning strategies, DNA fragments are initially generated and cloned in a non-specific manner, so that screening for the desired clone is carried out at the end of the process. Conversely, when specific DNA fragments are obtained by PCR or direct chemical synthesis, there is no need for subsequent screening.

WAYS OF GENERATING DNA FRAGMENTS




1. Non-specific generation of truly random fragments (by mechanical shearing or digestion with non-specific nucleases)
2. Through reverse transcription of mRNA into DNA
3. Highly specific amplification of a chosen piece of DNA by PCR
4. The use of synthetic DNA
5. Restriction endonucleases digestion

Dr. Sarob Khan

ENZYMES IN CLONING

Cutting and joining DNA



Dr.Sarookhani

RESTRICTION ENDONUCLEASES




- Enzymes
- Recognizes a short, symmetrical DNA sequence
- Hydrolyzes/cuts the DNA backbone in each strand
 - Specific site within that sequence
 - Foreign DNA is degraded into short fragments
- Finger printing

Dr.Sarookhani

RESTRICTION ENDONUCLEASES

2. Part of the restriction-modification defense mechanism against foreign DNA
3. Basic tools of gene cloning



Dr.Sarookhani

RESTRICTION ENDONUCLEASES

3 types

Type I

Type II

Type III



Dr.Sarookhani

TYPE II RES

- commonly used in cloning
- recognize and cut within (or immediately adjacent to) specific target sequences
 - generate specific fragments
- a small number
 - cut the DNA at a defined distance (usually only a few bases) away from the recognition site
 - limited applications
- requirement: Mg^{2+}



Dr.Sarookhani

CHARACTERIZATION AND IDENTIFICATION

1. The name of the organism from which they are obtained
2. Write in *italics*
 - The first letter of the genus
 - The first two letters of the species name
3. A suffix indicating the specific enzyme from that species

Dr.Sarookhani

CHARACTERIZATION AND IDENTIFICATION

- Example:
 - *Pst*I from *Providencia stuartii*
 - *Hae*I, *Hae*II and *Hae*III, three different enzymes, with different specificities from *Haemophilus aegyptius*

Dr.Sarookhani


THE PRODUCT OF REs DIGESTION

1. Products with protruding ends known as cohesive or 'sticky' ends
 - Fragments with unpaired single-stranded sequences either at the 5' or 3' ends



Dr.Sarookhani

THE PRODUCT OF RES DIGESTION



Dr.Sarookhani


THE PRODUCT OF REs DIGESTION

2. Products with blunt ends



Dr.Sarookhani

- Example




- Advantage: they can be joined to any other blunt-ended fragment
- Disadvantage: less efficiently ligated


Dr.Sarookhani

Enzymes	Recognition site	Number of bases	Ends generated	Source
<i>EcoRI</i>	G/AATTC	6	5' sticky	<i>Escherichia coli</i> RY13
<i>BamHI</i>	G/GATCC	6	5' sticky	<i>Bacillus amyloliquefaciens</i> H
<i>BglII</i>	A/GATCT	6	5' sticky	<i>Bacillus globigii</i>
<i>PstI</i>	CTGCA/G	6	3' sticky	<i>Providencia stuartii</i>
<i>XmaI</i>	C/CCGGG	6	5' sticky	<i>Xanthomonas malvacearum</i>
<i>SmaI</i>	CCC/GGG	6	blunt	<i>Serratia marcescens</i>
<i>Sau3A</i>	/GATC	4	5' sticky	<i>Staphylococcus aureus</i> 3A
<i>AluI</i>	AG/CT	4	blunt	<i>Arthrobacter luteus</i>
<i>NotI</i>	GC/GGCCGC	8	5' sticky	<i>Nocardia otitidis-caviarum</i>
<i>PacI</i>	TTAAT/TAA	8 Dr.Sarookhani	3' sticky	<i>Pseudomonas alcaligenes</i>

Restriction Enzyme Action of EcoRI



Dr.Sarookhani



OTHER ENZYMES IN CLONING

- Nucleases
- Ligases
- Phosphatase and Kinases
- DNA synthesizing enzymes

Dr.Sarookhani

NUCLEASES

Cutting/degrading DNA

ENDONUCLEASES

EXONUCLEASES

MULTIFUNCTIONAL NUCLEASES




Dr.Sarookhani

LIGATION

- The next stage in gene cloning
 - joining the DNA fragment to a vector molecule
 - DNA ligase




Dr.Sarookhani



PHOSPHATASES AND KINASES


- Removing or adding respectively phosphate groups
- Examples:
 - alkaline phosphatase: removes the 5' terminal phosphate from a DNA molecule leaving an OH group
 - polynucleotide kinase: adds a phosphate group to a free 5'-terminus
 - reverses the effect of alkaline phosphatase

Dr.Sarookhani




LINKERS AND ADAPTORS

- short synthetic DNA fragments that add new restriction sites to the end of a fragment




Dr.Sarookhani


L I N K E R S



A D A P T O R S



Dr.Sarookhani



TYPES OF CLONING VECTOR



Cloning limit
Purpose



Dr.Sarookhani

TYPES OF CLONING VECTOR (CLONING LIMIT)

Plasmid : 3-5 Kb

Phage : 8-23 Kb

Cosmids : up to 40 Kb

Bacterial Artificial Chromosomes : Up to 1 Mb

Yeast Artificial Chromosomes : Up to 1 Mb(eukaryot)

Retroviral : (for animal cells)

Phasmids (combination of PBR322 & Cosmid of
lambda phage)

Dr.Sarookhani

TYPES OF CLONING VECTOR (PURPOSE)

- Cloning vector (حامل های کلون)
 - reproduction of the DNA fragment (for gene study due to stability and easy)
- Expression vector
 - expression of certain gene in the DNA fragment

Dr.Sarookhani



CHARACTERISTICS OF A GOOD VECTOR

Stable

Easily detected

Easily isolated

Self-replicating


Multiple cloning sites

Small

Dr.Sarookhani

PLASMIDS


Dr.Sarookhani



Bacterial cells may contain extra-chromosomal DNA called plasmids.
Plasmids are usually represented by small, circular DNA.


(engineered plasmids can not transfer by itself)

Dr.Sarookhani




Schematic drawing of a bacterium with plasmids enclosed
(1) Chromosomal DNA (2) Plasmids.

Dr.Sarookhani



Selective marker is required for maintenance of plasmid in the cell. Because of the presence of the selective marker the plasmid becomes useful for the cell. Under the selective conditions, only cells that contain plasmids with selectable marker can survive. Commonly, genes that confer resistance to various antibiotics are used as selective markers in cloning vectors. For example, genes that render cells resistant to ampicillin, neomycin, or chloramphenicol are among commonly used selective markers.

Dr.Sarookhani



An example of a typical plasmid vector
(PUC19 has Amp resistant and Beta Gal. marker)

Dr.Sarookhani

Production of Whole Genomic library by phages or plasmids





Figure 10-23 Essential Cell Biology, 2/e. (© 2004 Garland Science)




Bacteriophage λ

(as a cloning vector)


Dr.Sarookhani

Bacteriophage lambda (I)




Transcriptional switches can
regulate cellular decisions

Dr.Sarookhani



Genes are clustered by function in the lambda genome



Dr.Sarookhani




Fig. 4.10 Map of the λ chromosome, showing the physical position of some genes on the full-length DNA of wild-type bacteriophage λ . Clusters of functionally related genes are indicated.

Cos site for circularization of genome within the E.coli

Dr.Sarookhani

Lysis or Lysogeny

Lysis: Infection by phage produces many progeny and breaks open (lyses) the host bacterium

Lysogeny: After infection, the phage DNA integrates into the host genome and resides there passively


No progeny


No lysis of the host

Can subsequently lyse (lyso~~geny~~)

Bacteriophage lambda can do **either**.


Dr.Sarookhani






Dr.Sarookhani


Vector construct




Dr.Sarookhani




In vitro packaging by cell lysate
Dr.Sarookhani





(Used for chromosomal walking)


Dr.Sarookhani



COSMIDS


- a plasmid containing a *cos* site
 - *cos* site: the sequence of bases of bacteriophage lambda that is cut asymmetrically during packaging, generating an unpaired sequence of 12 bases at each end of the phage DNA

Dr.Sarookhani



Structure of a cosmid

Dr.Sarookhani



YEAST ARTIFICIAL CHROMOSOMES (YACs) AND BACTERIAL ARTIFICIAL CHROMOSOMES (BACs)


Dr.Sarookhani

YACs AND BACs

- vectors capable of carrying very large cloned fragments
 - can accommodate genomic DNA fragments of more than 1 Mb (1 Mb=1000kb)
 - cloning the entire human genes
 - example: the cystic fibrosis gene ~ 250 kb
 - mapping the large-scale structure of large genome
 - example: the human genome project

Dr.Sarookhani

Vector transfer in yeast



Dr.Sarookhani

HOST

1. *E. coli*

well understood and conveniently manipulable organism

2. *Bacillus subtilis*

3. *Saccharomyces cerevisiae*

- frequently used eukaryotic host
- retains many of the convenient features of *E. coli*

4. *Picchia pastoris* (another yeast)

5. Cultured cells of *Spodoptera frugiperda* (an insect)

6- *Drosophila* for use with baculoviruses

7- cell lines

Dr.Sarookhani

Summary of E.coli Elite host Cells

Product	Genotype	Application
Elite Competent Cells (XL1-Blue)	recA1 endA1 gyrA96 thi-1 hsdR17 supE44 relA1 lac[F'proAB lac ^q ZΔM15 Tn10 (Tet ^r)]	Suitable for blue-white color screening and routine cloning applications.
Elite Competent Cells [BL21(DE3)]	F ⁻ ompT hsdS _B (r _B ⁻ m _B ⁻) dcm galλ(DE3)	Appropriate host for recombinant protein expression using T7-based expression vectors.
Elite Competent Cells (DH5α)	endA1 recA1 relA1 relA1 gyrA96 hsdR17(r _K ⁻ ,m _K ⁺) phoA supE44 thi-1Δ(lacZYA-argF)U169Φ80 Δ(lacZ) M15F ⁻	Suitable for cloning with large plasmid and cDNA library construction, and also allowing blue-white colonies selection.
Elite Competent Cells (JM109)	e14 ⁻ (McrA ⁻)recA1 endA1 gyrA96 thi-1 hsdR17(r _K ⁻ m _K ⁺)supE44 relA1 Δ(lac-proAB) [F ['] traD36 proAB lac ^q ZΔM15]	Appropriate for blue-white color and robotic screening (fast growing strain, forming visible colonies within 8~10 hours).

Dr.Sarookhani


Introducing vectors into host cells

- Transformation
- transfection
- Electroporation
- Nucleic acid gun
- Protoplast fusion
- Micro injection of DNA


Dr.Sarookhani

Figure 5.15 Two physical methods for introducing DNA into cells.

(a) Microinjection



(b) Transformation with microprojectiles



Dr.Sarookhani


Transformation is a mechanism of genetic transfer between bacteria in which the donor DNA exists cell-free in the recipient bacterium's immediate environment. DNA can be naturally released into the environment when cells die and subsequently lyse. Experimentally, DNA containing genes of interest, usually within a plasmid, can be introduced into the environment in order to transform bacterial cells. The ability of a recipient bacterium to take up free DNA and become transformed is known as competence. Some strains of bacteria are naturally competent. In others, competence is a brief physiological state during the exponential growth phase; in these bacteria, Ca^{2+} ions enhance the level of competence.

Dr.Sarookhani


TRANSFORMATION

- the most common hosts for simple cloning experiment: strains of *E. coli*
- treatment with solutions containing Ca^{2+} ions (sometimes Rb^{+} and Mn^{2+})
 - susceptible to take up exogenous DNA:
competent cells

Dr.Sarookhani



ELECTEROPORATION



- the diagram shows an electrical circuit diagram for a simple electroporation device.
- a cell suspension, such as of plant protoplasts or bacteria, is placed in the cuvette.
- the capacitor is charged by closing the right-hand switch. When the capacitor has been charged, the direct current pulse is discharged in the cuvette suspension by closing the left-hand switch.
- the DC pulse is thought both to disrupt temporarily the membrane and to electrophorese DNA into cells.

Dr.Sarookhani

TRANSFECTION

- equivalent to transformation

Dr.Sarookhani

SELECTION OF TRANSFORMANTS/MUTANTS

- sensitivity and resistance to chemicals
- requirement for certain compounds for growth (nutrient requirements)
- ability to use (breakdown) compounds
- plaque formation/type
- specific (by nucleic acid hybridization ((southern blotting)) or antibodies (western blotting)))

Dr.Sarookhani

INSERTIONAL INACTIVATION


Dr.Sarookhani

غیر فعال کردن ژن مقاومت به تتراسیکلین پلاسمید pBR322


- جایگاه آنزیم BamHI روی سکانس ژن مقاومت به تتراسیکلین پلاسمید pBR322 قرار دارد چنانچه DNA خارجی در جایگاه این آنزیم کلون شود و پلاسمید نو ترکیب در باکتری ترانسفرم گردد ، باکتری نسبت به تتراسیکلین حساس میشود. جایگاه آنزیم Pst I روی سکانس ژن مقاومت به آمپی سیلین قرار دارد.

Dr.Sarookhani


INSERTIONAL INACTIVATION



Dr. Sarookhan



REPLICA PLATING



Dr.Sarookhani




SELECTION OF TRANSFORMANTS IN BACTERIOPHAGE

Dr.Sarookhani


PLAQUE FORMATION/TYPE

- bacterial lysis and phage release occurs: a clear area is formed (clear plaque)
- lysogeny: turbid plaques

Dr.Sarookhani



3. Detecting the DNA sequence of a cloned gene with a probe (DNA hybridization)



Dr.Sarookhani






Fig. 6.13 Benton and Davis' plaque-lift procedure.

Western blotting



Screening can involve:

1. Phenotypic screening-
the protein encoded
by the gene changes
the colour of the
colony
2. Using antibodies that
recognize the protein
produced by a
particular gene



Dr.Sarookhani





Fig. 6.12 Granstein-Hogness method for detection of recombinant clones by colony hybridization.

Dr.Sarookhani


INSERTIONAL INACTIVATION OF LAC Z GENE (in Lac₋ Neg Bacteria)





سیستم Lac selection (Lac Operon)

- سلولي که حاوي وکتور کلونینگ با يك ژن LacZ ` بکر (دست نخورده) است میتواند بتا گالاکتوزیداز را سنتز کند . این آنزیم همراه با IPTG (Isopropyl-thiogalactoside) (که يك القا کننده است) X-gal (برموکلروایندولیل-گالاکتوپیرانوزید) را به ترکیب آبی رنگ برومو کلرو ایندولیل تبدیل میکند . این سلولها کلني آبی رنگ (فقط حاوي پلاسمید هستند) یا پلاک آبی رنگ (حاوي فاژ) تولید میکنند. سلولهاي نو ترکیبي که ژن LacZ ` آنها توسط Inserted DNA تخریب شده است نمیتوانند بتا گالاکتوزیداز را سنتز کنند و کلني هاي سفید یا پلاک هاي شفاف تولید میکنند

Dr.Sarookhani









- 1-prepro insulin method
- 2- separate A & B fragment production

Recombinant INSULIN production

(a) The artificial genes



(b) Synthesis of insulin protein

