



Anesthesia Machines and Systems

Circle Anesthesia Machine

- should provide for elimination of CO₂
- maintains directional gas flow
- allows artificial ventilation



Oxygen cylinder



“E” Cylinder is the size used with anesthesia machines

Compressed Gases (Oxygen)



- color coded
- 2000 psi (oxygen)
- 700 L. when full
- have pin index safety system
- “E” cylinder
- some gases support combustion (oxygen, nitrous oxide)


Nitrous Oxide




- compressed liquid
- 750 psi when liquid present
- analgesic
- relaxant
- more potent in humans
- used in 2-1 ratio with oxygen




Pressure gauge



Full oxygen tank = 2000 psi




Pressure reduction valve



Oxygen pressure reduced to 50 psi


Gas Pressure Regulators - Pressure Reduction Valves

- reduces high pressure in cylinder to 50 psi
- central supply oxygen pressure already reduced




Oxygen flowmeter

Measured in liters per minute




Flowmeter

- measures gas flow
- oxygen flows should meet or exceed metabolic requirement
- 2-3 ml/lb/min = metabolic requirement




Vaporizer (out of the circle)



Used with halothane, isoflurane, and sevoflurane

Vaporizer (in the circle)



Used with methoxyflurane (currently not in US)

Factors affecting output of vaporizers

- Barometric pressure
- Temperature
- Carrier gas flow rate
- Back pressure



Classification of vaporizers

- Precision
 - “Tecs”
 - “Matics”
 - Vapor
 - Saturation
- Method of output regulation
 - Variable bypass
 - Measured flow




Classification of Vaporizers


- Method of vaporization
- Location
- Temperature compensation




One way valves



A Tidal Volume is
estimated to be
■ 5 ml/lb




Rebreathing Bag



Recommended size = 30 ml/lb
or 6 x tidal volume

CO₂ absorption Canister



Recommended size = 2 X tidal volume

CO₂ absorption


- two formulations - soda lime and baralyme
- sodalime
 - 94% calcium hydroxide, 5% sodium hydroxide and 1% potassium hydroxide
- baralyme
 - 80% calcium hydroxide and 20% barium hydroxide

CO₂ absorption

- change absorbent after color changes 1/2 to 2/3
- change if granules are brittle
- “channeling” may be a problem



Pop off valve




Usually in the “open” position


Pop Off Valve

- can be either open or closed in a “closed” system (O₂ flow rate = 2-3 ml/lb/min)
- should be open in a “semi-closed” system
- should be either closed or partially closed if you are artificially breathing for the patient



Scavenger






- ### Bains Anesthesia System
- use oxygen flow rate of 100 ml/lb/min
 - used with small patients
 - no dead space
 - no resistance to breathing
- A vertical column of colored squares (purple, pink, blue, cyan, green, yellow, orange, red) is located on the right side of the list.


Artificial Ventilation (IPPV)

- ventilation rate = 8-20 breaths/minute
- tidal volume = 5-10 ml/lb
- inspiratory pressure = 10-30 cm/H₂O
- inspiratory time = 1 sec (2-3 sec in horses and cows)




Anesthesia Machines (Circle)

- Increased resistance
- dead space
 - the rebreathing tubes are not part of the dead space
 - length of E-tube outside of mouth is considered dead space
- oxygen flow rate
 - greater than 2-3 ml/lb/min



Anesthesia Machine (Closed)

- pop-off may or may not be closed
- flow rate is what determines classification
 - closed system = flow rate of 2-3 ml/lb/min
- economical
- no pollution
- heat conservation



anesthesia machine “closed”

- requires closer monitoring
- difficult to reanesthetize a patient that wakes up
- N₂O difficult to use
- buildup of toxic materials

