

Anesthesia Machine

The Anesthesia Machine

- ◆ You are the master of the machine
- ◆ You are responsible for checking the machine prior to each case

The Anesthesia Machine (con't)

- ◆ The primary cause of machine malfunction is failure to check
- ◆ Never start without the American Express items

What is the function of the anesthesia machine?

Functions of the Machine

- ◆ Convert supply gases from high pressure to low pressure
- ◆ Convert liquid agent to gas
- ◆ Deliver in a controlled manner

Functions (con't)

- ◆ Provide positive pressure for ventilation
- ◆ Alert the provider to malfunction
- ◆ Prevent delivery of a hypoxic mixture

Components of the Machine

- ◆ **Source gases**
- ◆ **Vaporizers**
- ◆ **Circuit**
- ◆ **Ventilator**
- ◆ **Scavenging system**

Safety Standards

- ◆ **1979 -- Standards set for all machines sold in the U.S.**
- ◆ **ANSI -- (American National Standards Institute)**
 - **Released 1979 standards**

Safety Standards (con't)

- ◆ **ASTM -- (American Society for Testing and Materials)**
 - Upgraded standards in 1988

The Generic Machine

- ◆ **2 sources of gas**
 - Pipeline 50 psig
 - Tanks
 - » Oxygen: 2200 psig
 - » Nitrous oxide: 745 psig
 - » Both reduced to 45 psig upon entering the machine

The Generic Machine (con't)

- ◆ Fail safe system (OFPD)
 - Stops flow if O₂ supply is lost
- ◆ Oxygen supply pressure alarm
- ◆ Second stage regulators
 - Reduces pressure to 14 psig

The Generic Machine (con't)

- ◆ Flow control valves
 - Regulate gas flow
 - Separates high and low pressure circuits
- ◆ Common manifold

The Generic Machine (con't)

- ◆ Vaporizer
- ◆ Outlet check valve
- ◆ Oxygen flush valve

Gas Sources

- ◆ Oxygen analysis is always required
- ◆ Pipeline
 - Enter at 50 psig
 - Gauge is on source side
 - DISS (Diameter Index Safety System)
 - » prevents gas swap

Gas Sources (con't)

◆ Side tanks

- Usually E cylinders
 - » Know pressure and volumes
- Enter at 45 psig
- Should be off unless in emergency use
 - » Prevents silent emptying

Gas Sources (con't)

◆ Pin index safety system

- Prevents tank swaps

Pin positions

Air 1-5

Oxygen 2-5

Nitrous oxide 3-5

Gas Sources (con't)

- ◆ Machine will use pipeline gas unless supply pressure drops below 45 psig

Fail Safe Devices

- ◆ Required by standards
- ◆ Stop flow of other gases if oxygen flow is interrupted
- ◆ Types
 - Threshold
 - Proportioning

Proportioning Systems

- ◆ Prevent delivery of less than 25% oxygen
- ◆ Either mechanical or pneumatic interface

Ohmeda Link-25 Proportion System

- ◆ Chain connects O₂ and N₂O flow control valves
- ◆ As N₂O is increased, the chain will turn O₂ control to maintain at least 25% O₂. Oxygen is increased

Ohmeda Link-25 Proportion System (con't)

- ◆ Maintains 3:1 ratio with combination of mechanical and pneumatic

Drager ORMC

- ◆ Pneumatic N₂O interlock
- ◆ Mobile shaft
- ◆ Slave control valve
- ◆ Pressure moves shaft and opens or closes slave valve

Drager ORMC (con't)

- ◆ N_2O flow is reduced to maintain 25% O_2
- ◆ Electrical contact provides alarm
 - Functional only in the O_2 / N_2O mode (not in the “all gases” mode)

Limitations of Proportioning Systems

- ◆ Wrong gas supply
- ◆ Defective operation
- ◆ Leaks downstream
- ◆ Inert gas administration

Flow Meter Assembly

- ◆ Controls and measures gas flow
- ◆ Thorpe tubes are tapered
- ◆ Indicator float is calibrated for specific tube
 - Density and viscosity differ
- ◆ Gas flows around float
 - Annular space

Flow Meter Standards

- ◆ Oxygen flow control knob
 - Physically different
 - Larger and projects further
 - Different shape
- ◆ All knobs are color coded
- ◆ Knobs are protected

Flow Meter Standards (con't)

- ◆ Low flow tubes for O₂ and N₂O
- ◆ Color coded flow tubes
- ◆ Thorpe tubes protected
- ◆ Tubes are not interchangeable
 - Float, tube and scale are single unit

Flow Meter Standards (con't)

- ◆ Note: Flow meters are located downstream from all safety devices except the oxygen analyzer.

Leaks

- ◆ Cracked tubes
- ◆ Faulty connections
- ◆ May create hypoxic mixture
- ◆ Oxygen is always downstream from other gases

Vaporizers

- ◆ Convert liquid anesthetic into a volatile inhalation agent
- ◆ Based on laws of physics
- ◆ You must memorize the chemical properties of the volatile agents

Applied Physics

- ◆ **Vapor pressure**
 - Dalton's law
 - Based on characteristics of agent
 - Varies with temperature

Applied Physics (con't)

- ◆ **Boiling point**
 - Vapor pressure equals atmospheric pressure
- ◆ **Latent heat of vaporization**
 - Heat required to change liquid into a vapor
 - Comes from liquid and environment

Types of Vaporizers

◆ Historic

- Copper kettle
- Vernitrol

◆ Modern

- Ohmeda Tec 4
- Drager Vapor 19.1

Ohmeda and Drager Characteristics

- ◆ Variable bypass
- ◆ Flow over
- ◆ Temperature compensated
- ◆ Agent specific
- ◆ Out of circuit

Copper Kettle and Vernitrol

- ◆ Measured flow
- ◆ Bubble through
- ◆ Non temperature compensated
- ◆ Multiple agent
- ◆ Out of circuit

Basic Design

- ◆ Gas enters vaporizer
- ◆ Flow is split
 - Majority is bypassed
 - Some enters vaporizing chamber
- ◆ Saturated gas leaves chamber
- ◆ Diluted by bypass gas
- ◆ Delivered to patient

Factors that Effect Output

◆ Flow rate

- Accurate at most flows
- Lower than dial setting at both extremes of flow

◆ Temperature

- Vapor pressure varies with temp
- Accurate at 20 - 35° C

Factors Effecting Output (con't)

◆ Intermittent back pressure

- Retrograde flow
- Higher than dial setting
 - » especially at low flows and high ventilator pressures

◆ Carrier gas composition

- N₂O causes transient drop

Vaporizer Interlock System

- ◆ Only 1 vaporizer can be turned on
- ◆ Gas enters only the “on” vaporizer
- ◆ Leak of trace gas is minimized
- ◆ Vaporizers are locked into the circuit

Vapor Pressures:

Isoflurane - 238

Enflurane - 175

Halothane - 241

Desflurane

- ◆ Requires special vaporizer
 - Vapor pressure 664
 - Pressurized, heated chamber
 - » 1550 mm / Hg prevents boiling

Vaporizer Hazards

- ◆ Misfilling
- ◆ Tipping
- ◆ Dual vaporizers on
- ◆ Leaks
- ◆ Free standing vaporizers

Misfilling

- ◆ Vaporizers are calibrated according to the vapor pressure of the agent
- ◆ If you fill with an agent with a higher v.p. -- overdose
- ◆ If you fill with an agent with a lower v.p. -- underdose

Anesthesia Circuits

Anesthesia Circuits

- ◆ **Link machine to patient**
- ◆ **Eliminate carbon dioxide**
- ◆ **Mapleson classification**
 - Many circuits in use
 - Modified Mapleson still in use
 - Know the current applications of modified Mapleson circuits

Types of Circuits

- ◆ **Basic circle system**
- ◆ **Mapleson Classification**

Basic components needed for delivery of anesthetic gases

Delivery Systems

- ◆ **Connection to patient**
- ◆ **Breathing tubing**
- ◆ **Unidirectional valves**
- ◆ **Breathing bag**

Delivery Systems (Cont'd)

- ◆ Pop-off valve
- ◆ Carbon dioxide absorption
- ◆ Bacterial filter

Circle System

- ◆ Allows rebreathing of anesthetic gases
 - lower FGF rates
 - Less pollution
- ◆ Requires CO₂ absorption
- ◆ Conserves heat and humidity

Advantages of Circle System

- ◆ Highly efficient
- ◆ Minimal dead space
- ◆ Conserves heat and moisture
- ◆ Minimal pollution
- ◆ Disadvantage - many places to leak

Components of the Circle System

- ◆ Fresh gas source
- ◆ Unidirectional valves
- ◆ Inspiratory & expiratory tubing
- ◆ Y-piece connector

Circle System Components (Cont'd)

- ◆ APL valve
- ◆ Reservoir bag
- ◆ CO₂ absorber

Rules for Circle System

- ◆ Unidirectional valve must be between patient & bag on both sides
- ◆ FGF cannot enter between patient & expiratory valve

Rules for Circle System (Cont'd)

- ◆ APL cannot be located
between patient & inspiratory
valve

Variations of the Circle System

Four Basic Circuits

- ◆ Open
- ◆ Semi-open
- ◆ Semi-closed
- ◆ Closed

Open Systems

- ◆ Insufflation
 - blow anesthetic gas over face
 - no direct contact
 - no rebreathing of gases
 - ventilation cannot be controlled
 - unknown amount delivered

Open Systems

- ◆ Open drop anesthesia
 - gauze covered wire mask
 - anesthesia dripped
 - inhaled air passes through gauze & picks up anesthetic

Open Systems (Cont'd)

- ◆ Open drop anesthesia (cont'd)
 - concentration varies
 - re-breathing may occur
 - environmental pollution

Semi-open Systems

- ◆ Breathing system which entrains room air
- ◆ Self inflating resuscitator system

Semi-closed System

- ◆ Gas enters from machine
 - part leaves via scavenger
- ◆ Circle system
- ◆ Bain system

Closed System

- ◆ Only enough gas enters to meet metabolic needs
- ◆ Scavenger is closed
- ◆ Closed circle system
- ◆ To-and-fro system

Closed System Anesthesia

- ◆ Technique not commonly used
- ◆ APL is closed and only enough O₂ is added to meet metabolic needs
- ◆ Anesthetic added based on square root of time
- ◆ Conserves anesthetic gas and eliminates pollution

The Scavenger System

- ◆ Releases excess pressure from the system
- ◆ Prevents operating room pollution
- ◆ Gases leave through APL
- ◆ May put too much negative pressure on the system

Systems Overview

Open System

- ◆ No reservoir
- ◆ No rebreathing

Semi-open System

- ◆ Has reservoir
- ◆ No rebreathing

Semi-closed System

- ◆ Has reservoir
- ◆ partial rebreathing

Closed System

- ◆ Has reservoir
- ◆ Complete rebreathing

Mapleson Breathing Circuits

- ◆ Early pioneers developed their own delivery systems
- ◆ Mapleson classified types of breathing devices

Mapleson Breathing Circuits (Cont'd)

- ◆ Mapleson circuits fall into which type of system?
- ◆ See Morgan p. 26, Table 3-1

Mapleson A

- ◆ FGI near bag
- ◆ Breathing tubing
- ◆ Expiratory valve near mask
- ◆ Volume of breathing tube should be as great as the tidal volume

Mapleson A

- ◆ Spontaneous ventilation
- ◆ High FGF flushes tubing between breaths

Mapleson A (Cont'd)

- ◆ Using “pop-off” enables controlled ventilation but also causes CO₂ rebreathing
- ◆ Current use?

Mapleson B

- ◆ Similar to A with FGI near expiratory valve
- ◆ System fills with FGF
 - inhaled by patient

Mapleson B (Cont'd)

- ◆ Exhaled gas forced out through expiratory valve
- ◆ Current use?

Mapleson C

- ◆ Similar to Mapleson B
- ◆ Shorter breathing tubing
 - less dead space
- ◆ Current use?

Mapleson D

- ◆ Long breathing tube
- ◆ FGI near mask
- ◆ Exhalation valve at distal end of breathing tubing
- ◆ Current use?

Bain Breathing Circuit

- ◆ Modified Mapleson D
- ◆ Tube within a tube
 - FGF tube within larger tube
- ◆ Mounts on anesthesia machine
- ◆ APL valve
- ◆ Connects to scavenger

Bain System

◆ Advantages

- compact, easy to handle
- warming of inspired gases
- partial rebreathing improves humidification
- APL controls system pressure
- ability of scavenging

Bain System Flow Rates

◆ Spontaneous ventilation

- 200-300 ml/kg/min

◆ Controlled ventilation

- | | |
|-----------------|--------------|
| – infants <10kg | 2 l/m |
| – 10 - 50 kg | 3.5 l/m |
| – > 60 kg | 70 ml/kg/min |

Bain System

- ◆ Depends on fresh gas flow to flush out CO₂
- ◆ Spontaneous ventilation
200 - 300 ml / kg / min
- ◆ Controlled ventilation
70 ml / kg / min

Mapleson E

- ◆ Exhalation tube is reservoir
–no bag
- ◆ FGI near mask
- ◆ Current use?

Mapleson F

- ◆ FGI near mask
- ◆ Breathing tubing/bag
- ◆ Expiratory valve at end of bag
- ◆ Current use?

Need To Know:

- ◆ Basic components
- ◆ Letters and names of systems currently in use
- ◆ Bain system
 - flow rates

Carbon Dioxide Absorption

- ◆ Allows rebreathing of anesthetic gases
- ◆ Review formulas from Chem / Physics
 - Know for Board exam

CO₂ Absorption (con't)

- ◆ Soda lime
 - 94% calcium hydroxide
 - 5% sodium hydroxide
 - 1% potassium hydroxide
 - silica to harden granules
 - ethyl violet as an indicator

CO₂ Absorption (con't)

◆ Baralime

- 80% calcium hydroxide
- 20% barium hydroxide
- ethyl violet as an indicator

CO₂ Absorption (con't)

- ◆ pH is extremely high
- ◆ Granule size
 - 4 8 mesh
- ◆ Water is required for chemical reactions to occur

CO₂ Absorber Incompatibility

- ◆ **Trichlorethylene**
 - dichloroacetylene
 - » neurotoxin
 - Phosgene
 - » pulmonary irritant
- ◆ **Sevoflurane**
 - » degrades in absorber

Ventilators Classified by:

- ◆ **Power source**
 - pneumatic
 - electric
 - both
- ◆ **Drive mechanism**
 - double circuit
 - driven by oxygen

Ventilator Classification (con't)

- ◆ **Cycling mechanism**
 - time cycled
 - pressure cycled
- ◆ **Bellows classification**
 - ascending / descending
 - » related to expiratory phase
 - Ascending is safer

Specific Ventilators

- ◆ **Review reading assignment**
- ◆ **Do not memorize technical data**
- ◆ **Note similarities and differences**

Ventilator Problems

- ◆ **Circuit disconnect**
 - Redundant alarms in place
 - Check APL valve
- ◆ **Occlusion**
- ◆ **Barotrauma**

Ventilator Problems (con't)

- ◆ **Leak in bellows assembly**
- ◆ **Mechanical problems**
- ◆ **Electrical problems**

Setting the Ventilator **(Things your mama didn't tell you)**

**Based on the principle that
PaCO₂ is directly proportional
to alveolar ventilation**

$$\text{AV} \times \text{CO}_2 = \text{AV} \times \text{CO}_2$$

(what you have) (what you want)

AV = alveolar ventilation

CO₂ = carbon dioxide

If you know 3, you can solve for the 4th

Patient weighs 150 lbs

R	10	20
TV	1000	500
MV	10,000	10,000
CO ₂	40	??

Alveolar Ventilation

- ◆ Minute ventilation minus dead space
- ◆ Dead space = 1 cc / lb

Ventilator Settings

- ◆ If rate is constant, then dead space is constant
- ◆ If you do not change the rate,
 $V_t \times CO_2 = X \times CO_2$

-
- ◆ You have $R = 8$, $V_t = 650$, $ETCO_2 = 40$. You want $ETCO_2 = 33$ and decide to leave the rate at 8. What new V_t is required to lower the $ETCO_2$ to 33?

$$V_t \times CO_2 = V_t \times CO_2$$

$$650 \times 40 = ?? \times 33$$

$$\text{New TV} = 788$$

Round off to 800 cc

Important concept

- ◆ PaCO₂ is directly proportional to alveolar ventilation
- ◆ If dead space is constant, alveolar ventilation is directly proportional to tidal volume.

Humidification

- ◆ Which takes more energy?
 - Humidification of dry gas
 - Heating cold gas

Humidifying a dry gas takes more energy than heating cold gas.

The Artificial Nose (Humidity Trap)

- ◆ Provides external heat and humidity
- ◆ More effective

Heated Humidifier

- ◆ More dangerous
 - Larger circuit volume
 - Increased circuit compliance
 - Thermal injuries

The Anesthesia Machine Check

- ◆ Required standard of care
- ◆ You are responsible for the function of your machine
- ◆ Follow the checklist

Machine Check (con't)

- ◆ Document “machine checked”
- ◆ Don't cut corners
 - Full check to start each day
 - Abbreviated check between cases

American Express Items

(Don't leave home without them)

Oxygen
Positive Pressure
Suction