

In The Name Of God

Site-specific PEGylation of an anti-CEA/CD3 bispecific antibody improves its antitumor efficacy

Presented by : fahimeh sadat mousavi alborzi

M.Sc student of medical biotechnology school of paramedical
sciences Qazvin university of medical sciences

Under supervision of : Dr. Farasat

Journal information

Journal impact factor:4.370 (5 year impact 5.14)

Indexed in :

PubMed and PubMed central(Int J Nanomedicine)
MedLine


American chemical society's chemical abstracts
service(CAS)

Science citation **Index** Expanded &...

Article information:

Author :Haitao and et al

Published:2018 may 29



Index

Introduction.....	4
Material and method.....	8
Result.....	24
Discussion.....	32
Conclusion.....	35

Introduction

Antibodies have attracted increasing interest as **therapeutics** due to their unique **biophysical** characteristics, such as **good solubility**, **stability**, **high specificity**, and **prolonged plasma concentration**.

bispecific antibodies have been developed as a powerful approach to cancer immunotherapy by engaging **immune cells** to target **cancer cells**.

Bispecific antibodies have two different antigen-binding sites:

- 1) recognizes tumor cells
- 2) recognizes immune cells, usually T cells or natural killer(NK) cells.

Nb Derived from natural camel heavy-chain only antibodies
single-domain antibodies lack the first constant (CH1) domain and light chain and are consequently referred to as nanobodies or VHHs

Nanobodies are small sized

more stable than c ScFvs or BiTEs (bispecific t-cell engager)
a good scaffold for constructing bispecific antibodies

- S-fab was developed -can bind to diverse epitope-be synthesized in prokaryotic expression systems.
- S-Fab has potent in vitro **cytotoxicity**
- suppresses **cancer progression** in in vivo studies
- have a short plasma **half-life**

To improve the in vivo half-life of proteins:

- 1) conjugation to **polyethyleneglycol** (PEG; PEGylation)
- 2) N-(2-hydroxypropyl) methacryl-amide (**HPMA**) copolymers
- 3) proteins (such as **albumin**)
- 4) **polyglutamic acid** and **PASylation**

Other benefits of PEGylation

reducing the antigenicity and immunogenicity of the conjugated proteins
improving protein solubility
enhancing the proteolytic resistance of therapeutic proteins
decreasing toxicity, and improving the thermal and mechanical stabilities

In this study, we explored **thiol site-specific PEGylation** to improve the **half-life ($t_{1/2}$)** of the CEA-S-Fab bispecific antibody, which is an antiCEA/CD3 bispecific S-Fab

A functionalized **20** kDa linear PEG was conjugated to S-Fab.

To retain the properties of CEA-S-Fab, the site for conjugation was designed to be distant to both **CEA- and CD3-binding sites**.

Material and methods



The 20 kDa linear methoxy PEG maleimide (MAL-PEG-OMe; Lot A3115)

Human colon adenocarcinoma cell line **LS174T** and human breast cancer cell line **SKOV3**
(NOD–SCID) mice (18–22 g)

Adult male Sprague Dawley (SD) rats(250–300 g) that were used for PK analysis

Methods

1) S-Fab design and protein purification



To produce S-Fab, two plasmids encoding the **VH-CH1-VHH** and **VL-CL** polypeptides were co-transformed into **BL21** (DE3, codon plus) competent cells with the appropriate antibiotics.

When the absorbance (**OD600**) of the cell culture reached 0.8
0.2 mM isopropyl- β -d-thiogalactoside (IPTG) was added to induce protein expression.

Cells were cultured at
16°C for another 40 h before harvesting.

periplasmic extraction was performed by re-suspending the cell pellets in a pre-cooled sucrose solution (20 mM Tris-HCl, pH 7.5; 25% [w/v] sucrose; 1 mM EDTA).

After **incubation on ice**, the suspension was **centrifuged** and the supernatant fraction was **collected** as the **sucrose** fraction.

The pellet was re-suspended in a chilled periplasmic solution (5 mM MgCl₂) and centrifuged .
supernatant was gathered as the periplasmic fraction.

The S-Fab was purified from the combined sucrose and periplasmic fractions by a two-step :

- 1) by immobilized Ni-NTA affinity chromatography
- 2) by an IgG-CH1 affinity matrix

Gel filtration analysis was performed using a Bio-Rad FPLP system

analysis under reducing conditions.

The resulting fractions were visualized by **Coomassie blue** staining.

Protein markers (Lot MWGF200; Sigma-Aldrich Co.) were loaded as standard controls for gel filtration analysis.

2)Conjugation of S-Fab to PEG

- S-fab in 5ml **pbs** & 1 mM **tris 2-carboxy ethyl phosphine** were mixed and incubated to obtain reduced s-fab fragment
- To explore the optimal molar ratio of **MAL-PEG-Ome** and **S-Fab** in the PEGylation
- a series of reactions with the molar equivalents of PEG:S-Fab of **0:1, 10:1, 20:1, 40:1, and 60:1**.
MAL-PE - OMe was dis-solved in sterile water to obtain a working concentration of 20 mg/mL (1 mM).
- **PEGylation of S-Fab was carried out** by mixing MAL-PEG- OMe (at the working concentration) with reduced S-Fab and shaking at 22°C for 2 h.

- Result ➡ samples ➡ 12% reduced ➡ **SDS-PAGE**
or non reduced
- ➡ **western blot** ➡ transfer polyvinylidene fluoride
bloking (skimmed milk) ➡ membrane were **incubated** with
antiflag HRP(1:2000 heavy chain) and anti-His IgG(1:3000 light-
chain) in 5% skim milk ➡ after **washing** with tris-buffered
saline & Tween 20 buffer ➡ secondary **Ab** (HRP-conjugate
IgG1:13000) Incubated with light-chain ➡ membrane were
developed pierce's west pico chemiluminescence

3) Purification of PEG-S-Fab using fast protein liquid chromatography

PEG-S-Fab was purified using an AKTA™ avant25 fast protein liquid chromatography purification system (GE Healthcare Bio-Sciences Corp.) and a Superdex 10/300 GL column at a flow rate of 0.8 mL/min. The column was first equilibrated with two column volumes (CVs) of distilled water and two CVs of PBS before applying the samples. All the collected fractions were analyzed by Coomassie blue and barium iodide complex staining after SDS-PAGE under reducing conditions. The fractions of the purified PEG-S-Fab were pooled together for further studies.

4) Human CD3 T-cell isolation

Human PBMCs were prepared from healthy donors using Ficoll gradient centrifugation as previously described.^{8,11} T cells were isolated from PBMCs using an EasySep™ Human CD3 Positive Selection Kit. Isolated T cells were cultured in complete Roswell Park Memorial Institute 1640 medium with 10% FBS and 1% penicillin/streptomycin at 37°C in a 5% CO₂ humidified incubator before cytotoxicity assays.

5)Flow cytometric analysis

1)Detect CEA → cell line **LS174T & SKOV3** → centrifugation

→ **washing** ice-cold PBS+0.2% BSA → 1×10^6 LS174T & SKOV3 or

5×10^5 TCELL → S-Fab & PEG_sfab **incubated** on **ice** for 1h →

washing twice with ice-cold PBS +BSA → Anti-CD3 fluoresceine

isothiocyanate (**OKT3**) as +control **for CD3+antigen-binding analysis** .

→ Goat antihuman IgG (H +L)-**AlexaFluor 488 antibody** added → cells
were **incubated** on **ice** for 1 h → After **washing** the cells twice,. → **Flow**
cytometric

6) Immunofluorescence assay

- 1) LS174T & SKOV3 were **plated** 30 mm confocal glass bottom dishes
- 2) **Washed** with cold PBS three times
- 3) Fixed cells were **incubated** with S-fab or peg-s fab + IgG (H+L)AlexaFluor488 antibody for 2 h at 4°C.
- 4) Cell nuclei were **counterstained** with 4',6-diamidino-2-phenyl indole (DAPI).
- 5) **washing** with PBS,
- 6) the samples were **examined** under a scanning confocal microscope
- 7) **analyzed** using the Olympus FV31S-SW_V2.1 software.

7) In vitro cytotoxicity

- 1) LS174 & SKOV3 → target cells
- 2) CD3+ T cell without prior stimulation → effector cells
- 3) Cytotoxicity assay → 96-well microplate -triplicate-seeding 5000 target cell in 100ML media
- 4) 6 h incubation an equal volume of CD3+ Tcell add to each well an E:T ratio of 10:1
- 5) add a series of concentrations (0.033, 0.1, 0.33, 1, 3.3, 10, 33, and 100 nM) of S-Fab or PEG-S-Fab
- 6) After a 72 h incubation, the cell viability was evaluated

8) PK assay to determine the in vivo half-life

1) Specific pathogen free male **SD rats** (250–300 g) were used for the PK assay.

2). **S-Fab** , **PEG-S-Fab**, or a volume equivalent of the **vehicle solution** **PBS** was administered through the **caudal vein**.

3) **0, 0.5, 1, 2, 4, 8, 16, 24, 36, 48, 72, 96**, and **144** h after administration
. ➡ Sample blood ➡ orbital vein ➡ using capillary

4) Blood sample ➡ centrifugation ➡ stored at -80

5) S-Fab & PEG-S-Fab in the plasma samples were quantified using **ELISA**.

6) 100 μL aliquot of 6D6 (mouse antihuman IgG Fab antibody; 1.0 $\mu\text{g}/\text{mL}$ in PBS) was used to coat each well of a 96-well ELISA microplate for 2 h at 37°C.

7) The wells were blocked with blocking buffer (PBST containing 1% BSA)

8) Samples and standards (100, 80, 50, 40, 30, 20, 10, 5, 1, and 0.1 $\mu\text{g}/\text{mL}$) in
→ blocking buffer → (S-Fab)/PBS 1:10 → important, to avoid matrix effects in the assay

9) 1:3 dilution plasma

10) Sample added in triplicate and incubated

11) each well → washed with PBST → added secondary

Ab(mouse monoclonal anti-flagM2-peroxidase [HRP] antibody at 1:500 dilutions) per well

12) Five time washed → 100 µL **tetramethylbenzidine** substrate solution was add to each well

13) Incubation

14) 100 µL of 2 M **H₂SO₄** Was added to stop the reaction

15) The absorbance was detected at **450 nm** using a Tecan ELISA microplate reader.

16) The serum elimination $t_{1/2}$ and clearance were calculated with **3P97 PK software**

.

9) Plasma stability analysis

S-Fab and PEG-S-Fab stabilities were assessed in human fresh plasma over 2 weeks. Briefly, S-Fab and PEG-S-Fab were diluted with human fresh plasma (without platelets), which generated an initial concentration of 100 $\mu\text{g/mL}$. Simultaneously, a vehicle control with only plasma was established. The samples were incubated at 37°C for 2 weeks. At the time intervals of 0, 24, 48, 72, 96, 168, 264, and 336 h, 40 μL samples were collected and stored directly at -80°C until further analysis. The samples were thawed on ice and centrifuged at 18,500 $\times g$ for 10 min at 4°C. The supernatant was subjected to electrophoresis in 12% reducing SDS-PAGE (5 μL of sample per well). After electrophoresis, a Western blotting assay was performed to analyze the protein levels.


10) In vivo **tumor** growth inhibition

- 1) using **NOD-SCID** mice engrafted subcutaneously with LS174 cell
- 2) Cells **harvested** and **washed** once with **PBS** and mixed with human **PBMCs** that were freshly isolated from healthy donors.
- 3) Mixtures of LS174T cells and human PBMCs were **subcutaneously** injected into the right flank of NOD–SCID mice
- 4) One hour after engraftment, S-Fab & PEG-S-Fab or the vehicle control (PBS) were injected **intraperitoneally**.
- 5) The tumor volume was measured with **calipers** in two perpendicular dimensions
- 6) All data were expressed as the mean \pm **SEM** for each group, and differences between groups were determined by a two-way analysis of variance using the **GraphPad Prism 5** software

Result



1) Bispecific S-fab construction and purification

- S-fab was constructed with the following specification
- Peptid(CGGGGC) added to C-terminal of anti-CD3 CL-VL
- Cloned into PET26b & PET21a
With signal peptide for pelB



- C)
 - Light chain
25KD
 - Heavy chain
40KD
- (under reducing condition 2-mecapto ethanol)

A major binding
130KD observed
(under
nonreducing)=
Dimerzation s-
fab



d)Gel filtration analysis
Intact s-fab Ab 130KD

2) Conjugation of s-fab to PEG




Figure 2 PEGylation of bispecific S-Fab using 20 kDa linear MAL-PEG-OMe.

Notes: S-Fab was reacted with the functionalized PEG at a series of ratios (0:1, 10:1, 20:1, 40:1, and 60:1 molar equivalents of PEG:S-Fab). After shaking at 22°C for 2 h, the samples were examined by reducing SDS-PAGE (5 μ L per sample). **(A)** Coomassie blue staining of S-Fab with different molar ratios of PEG:S-Fab during PEGylation. **(B)** Barium iodide complex dye staining for S-Fab with different molar ratios of PEG:S-Fab. **(C)** Western blotting assay using an anti-FLAG antibody to detect the VH-CH-CEA chain. **(D)** Western blotting assay using an anti-His6 antibody to detect the VL-CL chain. Numbers 0, 10, 20, 40, and 60 represent the molar ratio of PEG:Fab at 0:1, 10:1, 20:1, 40:1, and 60:1, respectively. M (kDa), molecular weight markers (kilodalton).

Abbreviations: CEA, carcinoembryonic antigen; MAL-PEG-OMe, methoxy PEG maleimide; PEG, polyethylene glycol; S-Fab, single-domain antibody-linked Fab; SDS-PAGE, sodium dodecyl sulfate polyacrylamide gel electrophoresis.

3) Purification of PEG-S-Fab using fast





Figure 3 Purification of PEG-S-Fab.


Notes: Gel filtration was used to fractionate the PEGylation mixture (fractions 1, 2, and 3). **(A)** Chromatogram of PEG-S-Fab. **(B)** Coomassie blue staining of the fractions 1, 2, and 3 from PEG-S-Fab purification. **(C)** Barium iodide complex staining of the fractions 1, 2, and 3 from PEG-S-Fab purification. M (kDa), molecular weight markers (kilodalton).

Abbreviations: PEG, polyethylene glycol; PEG-S-Fab, PEGylated S-Fab; S-Fab, single-domain antibody-linked Fab; UV, ultraviolet.


4) PEG-S-Fab can bind the tumor antigen CEA and CD3 T cells



6) PEGylation prolongs the in vivo half-life of s-fab



7) PEG-S-fab induces more potent in vivo antitumor activity



Discussion

- In this study, we demonstrated that the site-specific PEGylation can improve the therapeutic profiles of the bispecific
- S-Fab showed potent antitumor activities in vitro and in vivo
- As a Fab format, S-Fab offers advantages over full-length antibodies, including **better penetration** into tumor tissues and **simple production** of the antibody using inexpensive **prokaryotic expression systems**.
- the short in vivo half-life of S-Fab limits its potential as a therapeutic agent.

In the case of antibody fragments, PEGylation has been shown to **lengthen** their circulating **half-life**, enhance their proteolytic resistance of therapeutic proteins, and **reduce** their **immunogenicity**.

PEGylation potentially decrease enzymatic proteins or the binding affinities of Fabs.

Flow cytometry analysis showed that PEGylation had a slight effect on the binding of S-Fab to tumor cells and T cells and slightly reduced the in vitro cytotoxicity

PEGylation increased the circulation $t_{1/2}$ of S-Fab by 12-fold in vivo

. enhancement of the in vivo half was due to the **hydrodynamic radius** and **enhanced molecular size** of the PEGylated conjugate leading to **reduced renal clearance** and **increased resistance** of the PEGylated conjugate to proteolysis.

In mouse models, PEG-S-Fab showed **more potent antitumor activity** than S-Fab

Conclusion

this study **supports the use of PEGylation** to extend the circulating half-life of the bispecific antibody S-Fab and enhance **its therapeutic potential** in the clinic. This strategy can be potentially applied to a broad range of immunotherapies.

با تشکر از توجه شما

